



Speakers' setup

Step 1: Select four long female-to-female jumper wires and connect one end of each to the pins of the speakers.




Step 2: Remove the square female heads from the other ends of the wires.


Step 3: Solder the ends of two jumper wire (one from each speaker) together along with one small female-to-female jumper wire (after removing its square female head).



Step 4: Cut a piece of heat shrink tube long enough to cover the exposed wires.




Ready! Time for testing!



Testing the speakers


Step 5: Use the oscilloscope and connect one end (*measuring clip*) to one wire of the speaker and the grounding clip to the other wire.



Step 6: Connect your speakers to your Raspberry Pi Pinout. First, go to the `EthoPy/Interfaces/ RPPorts.py` and look at the "class `RPPorts(Interface)`". Here, you see which channel corresponds to the RP GPIO pins based on what you want to connect. So, for 'sound' is channel GPIO 13.

```
class RPPorts(Interface):
    channels = {'Odor': {1: 24, 2: 25},
               'Liquid': {1: 22, 2: 23},
               'Lick': {1: 17, 2: 27},
               'Proximity': {3: 9, 1: 5, 2: 6},
               'Sound': {1: 13},
               'Sync': {'in': 21, 'rec': 26, 'out': 16},
               'Opto': 19,
               'Status': 20}
```

Connect one wire of your speaker to **GPIO13** and the other wire to **ground**.



Step 7: Open your RP and run a task with auditory stimuli.